



Gantry Stand Selection Guide

Ultimate Flexibility in Positioning Micromanipulators

GANTRY STAND SELECTION GUIDE

TALL GANTRY

For use with an inverted microscope.


Independent tower solution when not using our MD manipulator or stage platforms.

STANDARD GANTRY

For use with an upright microscope mounted on a translator

SHORT GANTRY


For use with an upright microscope without a translator.



BASE PLATE OPTIONS


SQUARE BASE OPTIONS

Short: **505796**
Standard: **505996**
Tall: **505997**
6 in L x 6 in W x 0.75 in H




NARROW BASE OPTIONS

Short: **505990**
Standard: **505991**
Tall: **505992**
3 in x 5 in x 0.75 in



THIN BASE OPTIONS

Short: **505993**
Standard: **505994**
Tall: **505995**
5 in x 5 in x 0.3 in



BASEPLATE VERTICAL DIMENSIONS

Stand	Tall		Standard		Short	
SQUARE OR NARROW	11.71 in 29.74 cm	15.25 in 38.74 cm	8.71 in 22.12 cm	12.25 in 31.12 cm	6.71 in 17.04 cm	7.25 in 18.42 cm
THIN	11.26 in 28.60 cm	14.80 in 37.59 cm	8.26 in 20.98 cm	11.80 in 29.97 cm	6.26 in 15.90 cm	6.80 in 17.27 cm